

Stretch Film

StretchFlex® PCR

IPG®

TECHNICAL DATA SHEET

Test	Values (listed by film thickness)		
	Reference Methods	63 gauge 16.0 micron 0.63 mil	80 gauge 20.3 micron 0.80 mil
Ultimate Tensile MD (psi)	ASTM D882	6500	6200
Ultimate Tensile TD (psi)	ASTM D882	3500	3000
Elongation MD (%)	ASTM D882	540	600
Elongation TD (%)	ASTM D882	850	825
Load Retention - Average (3 top, 3 bottom, 50% overlap @ 75%)	Lantech Tool	10.00	12.00
Tear Time (seconds)	ESTL	10.5	12.0
Pre-Stretch (% at break)	IPG/Highlight	230	250
Cling (in to out, grams at 100% pre-stretch)	IPG method	100 grams	100 grams



100 Paramount Drive, Suite 300
Sarasota, FL 34232 | USA
Customer Service: 888.898.7834
www.itape.com | info@itape.com

While we believe them to be reliable, the statements and information herein are only for general guidance and are not warrants or guarantees for accuracy and completeness. The user must, by test or otherwise, determine suitability for this purpose. There is no warranty of fitness for a particular purpose. Our standard term and conditions of sale apply exclusively to all orders, and all liability for damages of any kind, including consequential, exceeding purchase price is excluded. No one is authorized by us to make oral warranties. We reserve the right to make changes without notice or obligation in our products and publications.

EFFECTIVE: 10/25

